



Soleoline

SOLEO LINE

Ultrasound, electro- and combination therapies

Soleo Sono, Soleo Galva und Soleo SonoStim

These three therapeutic devices from the Soleo-line series are innovative, powerful and very easy to use thanks to their small size. Ideal for mobile use. Whether in your practice, during house calls, or at sporting events, you can use them wherever your patients are.

Extremely easy to use

- A large, colour display with one-touch operation
- Easy programme selection based on treatment goal
- Easy to use thanks to intuitive menu navigation



Expertise for any medical problem

Clearly presented, recommended treatments for various therapies available at all times.

Tailored to your practice

- The programme's customised start settings save you precious treatment time.
- Change pre-set programmes and save them to fit your patients.
- Your custom settings are saved to an SD card.

The comprehensive design of the new *Soleoline* generation combines aesthetics with functionality, ergonomics and durability.

Simply beautiful.

A sleek case and an attractive, capacitive glass touchscreen ensure that *Soleo Sono*, *Soleo Galva* und *Soleo SonoStim* really stand out in your practice.

Beautifully simple.

The user interface of the new *Soleoline* generation has been redesigned to be clear and self-explanatory. It's even fun to use with every application.

Intelligent design –
sets new standards



Soleo Somo combines the classic possibilities of modern ultrasound therapy with the new SonoSwing process.

Helps wherever the pain is

The ergonomic, multiple frequency transducer can give off two frequencies: 0.8 MHz (for penetration) and 2.4 MHz (to treat the surface).

Unique

The patented SonoSwing process allows you to adjust the percentage of both frequencies. Due to the adjustable mix of both frequencies you can choose the depth of the optimal heat development.

Time-saving and safe

- Select treatment based on symptoms
- Recommended treatments for 50+ symptoms
- Up to 120 favourite and customised programmes can be stored

The new dimension of ultrasound therapy

SOLEO SOMO

For the full spectrum of modern electrotherapy

With a wide variety of current types, special sports and muscle rehabilitation programmes Soleo Galva offers a comprehensive range of treatments.

Two channels

The number of different varieties in the 2-channel mode increases your range of therapies.

Twinstim: 1 programme / 2 channels

Duostim: 2 programmes / 2 channels

Sports therapy

Special additional current types and programmes allow you to expand the different types of sports therapy you offer.

VacoS: Electrotherapy plus suction massage

VacoS allows you to apply electrodes with ease. Your patients are sure to enjoy the therapeutic effects of a suction massage during electrotherapy.



SOLEO SONO STIM

SOLEO SONO STIM

Ultrasound and electrotherapy all in one

Soleo *SonoStim* does everything that Soleo *Sono* and Soleo *Galva* do, and it also allows you to combine ultrasound therapy with electrotherapy. With its innovative technology, more than 300 indications, this mobile and high-tech device has a huge range of applications.



Simultaneous therapy

The simultaneous delivery of ultrasound and electrotherapy intensifies the therapeutic effects of each.

Ultrasound and electrotherapy combined
- Really easy to use

Due to an adjustable mix of both frequencies
you determine the optimal depth of heat Generation.

High-performance therapy

Technical Data



Ultrasound therapy version

Ultrasound head

Frequencies

Maximum power emission

User guide

Sono

1 socket for one 5 cm² multi-frequency watertight head, (1 cm² as an optional extra)

0.8 MHz: deep effect / 2.4 MHz: superficial effect

3 W/cm² for the 5 cm² head

More than 50 therapy recommendations

Electrostimulation version

Channels

User guide

Galva

2 independent channels (programmes, intensities, parameters)

More than 120 therapy recommendations

Combined version

General features

User guide

SonoStim = Sono + Galva

Same features as **Sono** and **Galva**

More than 300 therapy recommendations

Soleoline

SD memory card

Dimensions

For all versions

120 favourites entries

326 mm x 134 mm x 243 mm / 2,1 kg

Vaco equipment

Operating modes

Control

Dimensions

VacoS for Galva or SonoStim version

Continuous and pulsed

via the touch screen

322 mm / 234 mm / 130 cm / Weight 3 kg

Zimmer MedizinSysteme GmbH
Junkersstraße 9
89231 Neu-Ulm, Germany
Tel. +49 7 31. 97 61-291
Fax +49 7 31. 97 61-299
export@zimmer.de
www.zimmer.de

