



Micro*Pro*

MICRO PRO

The safe microwave therapy



With the *MicroPro* you have a powerful microwave therapy device at your disposal, the functional principle of which has proven itself in everyday practice for many years.

The microwave energy is transferred to the tissue with the help of a radiator and is then converted into heat. Depending on the type of treatment, the user can choose between different types of radiator: large-area, round, trough-shaped or point-type. A freely adjustable swivel arm is provided to focus the radiator on the body part to be treated.

To guarantee ease of use as well as convenient and efficient operation and adjustment of parameters, the unit features a 12"-display with capacitive touch-screen.

The treatment time can be set digitally between 1 and 30 minutes. At the end of the treatment, the set power is automatically switched off and the end of the treatment is signalled visually and acoustically. By means of the heat-effective output power value that is constantly displayed, the user can exactly adjust parameters and has full visual control over the therapy.

In addition to the continuous operating mode, the *MicroPro* offers a pulse function for the treatment of deeper tissue layers with simultaneous low thermal surface exposure. Thus *MicroPro* can be used across a broad range of applications. As a result, *MicroPro* has a wide range of applications in hospitals, among doctors and physiotherapists.



Additional features

- Regulated fans
- Modern user interface
- Modern housing
- Favourites can be saved
- Favourites can be exported via USB and SD-card

MICRO PRO

MICRO PRO



Technical data

MD

HF nominal power	200/200 W (CW/pulse)
Operating frequency	2.45 GHz \pm 1.5 %
Power indication	Effective power
Timer	1–30 min. (1 min steps)
Power supply	230 V, 50/60 Hz (Modifiable to 115 V)
Power consumption	900 VA
Protection class	I/B (EN 60601-1)
Equipment class	Ila
Dimensions	450 x 1110 x 530 mm (B x H x T)
Weight	44 kg

Radiators



Longitudinal-field radiator



Circular-field radiator



Cradle radiator



Focus radiator

Distributor:

Zimmer MedizinSysteme GmbH
Junkersstraße 9
89231 Neu-Ulm, Germany
Tel. +49 (0)7 31. 97 61-291
Fax +49 (0)7 31. 97 61-299
export@zimmer.de
www.zimmer.de

Zimmer
MedizinSysteme

